VETIGEL®

HAEMOSTATIC GEL

Rx ONLY, ANIMAL USE ONLY. NOT FOR HUMAN USE.

INSTRUCTIONS FOR USE

Each sterile, single-use 5-ml haemostatic gel syringe is individually pouched.

Caution: For animal use only. This device is restricted for use by or on the order of a veterinarian.

DESCRIPTION

VETIGEL® Haemostatic Gel ("VETIGEL") is a single-use, prescription haemostatic gel for veterinary use only.

VETIGEL® is intended as an adjunctive haemostatic device.

VETIGEL® is supplied as individually pouched 5-ml haemostatic gel syringes (sodium alginate and poly (N-acetyl-D-glucosamine, D-glucosamine) hydrogel).

The haemostatic gel is an opaque, tan, colloidal hydrogel.

The two polymer components of VETIGEL®, sodium alginate and poly (N-acetyl-D-glucosamine, D-glucosamine), are sterilised individually by irradiation, the haemostatic gel is manufactured and filled into syringes and the pouched filled syringe is terminally sterilised with gamma irradiation.

All components of VETIGEL® are non-animal derived.

Due to the viscous nature of the material, benign air bubbles may be present in the syringe (no impact on safety or efficacy).

VETIGEL®'s polymer components are naturally-derived and thus may differ in colour within any given syringe or over time. Particles may appear dark at times.

INTENDED USE/INDICATIONS

VETIGEL® Haemostatic Gel ("VETIGEL") is a veterinary haemostatic gel for animal use only. VETIGEL® is indicated as an adjunctive haemostatic device to assist when control of bleeding by pressure, ligation and other conventional procedures is ineffective or impractical.

CONTRAINDICATIONS

- The safety and effectiveness of VETIGEL® for use in ophthalmic procedures have not been established.
- VETIGEL® is not intended to be used as a wound-closure device.

WARNINGS

- Do not use VETIGEL® in *Equus* due to pyrogen sensitivity.
- Do not inject VETIGEL® intravascularly due to risk of embolisation.

PRECAUTIONS

- VETIGEL® can only be used by or under the direction of a veterinarian.
- Do not use VETIGEL® if the packaging seal is torn and/or damaged.
- VETIGEL® must not be used beyond its expiry date.
- VETIGEL® is manufactured and supplied as a one-time-use medical device. It cannot be sterilised before use or after its contents have been opened or exposed.
- VETIGEL® is not a substitute for good surgical technique.
- The use of VETIGEL® is not recommended in the presence of an active infection. If symptoms of infection or abscessation develop where VETIGEL® has been deployed, revision surgery may be necessary to remove the infected material and allow for drainage.
- Degradation profile for VETIGEL® Haemostatic Gel is > 30 days, definitive study is ongoing.
- VETIGEL® was assessed for its biocompatibility per ISO 10993 (human medical device biocompatibility standards), interpreted for veterinary use.

MECHANISM OF ACTION


VETIGEL® is comprised of a proprietary blend of polyanionic and polycationic polysaccharides. Sodium alginate is the primary polyanionic polymer, and poly (N-acetyl-D-glucosamine, D-glucosamine) is the primary polycationic polymer.

Sodium alginate forms a hydrogel in which poly (N-acetyl-D-glucosamine, D-glucosamine) particles are uniformly dispersed. When directly applied to a source of bleeding, the haemostatic gel rapidly adheres to the wound site. The haemostatic gel forms a mechanical barrier that stops the flow of bleeding and allows the body to create a natural clot without the need for applied pressure.

DIRECTIONS FOR USE


- Retrieve one VETIGEL® syringe from box. Though VETIGEL® does not require preparation, it is recommended that prior to beginning surgery, one VETIGEL® syringe is pre-emptively taken out of the box and kept readily accessible in case of need.
- Using aseptic technique, open pouch and deploy syringe into sterile field. When ready to use, unscrew syringe cap and prime syringe by slowly depressing plunger until gel reaches the nozzle. Wipe any excess gel that comes out from the nozzle with sterile gauze.

Prime the Syringe

- Identify the source of bleeding.
- Swab the wound/site of bleeding with sterile gauze. Remove as much blood as possible in order to allow the haemostatic gel to come into direct contact with the source of bleeding.


Remove Excess Blood

- Remove gauze and apply the full syringe contents at the deepest part of the lesion holding the syringe perpendicular to the wound surface. Once the gel has filled the wound cavity, start slowly withdrawing the syringe from the defect while continuing to expel the remainder of the gel. When the gel is flush with the surrounding tissue, shear the gel by bringing the nozzle of the syringe away from the defect and downward toward the adjacent tissue in a wiping motion to ensure the gel anchors to the tissue and does not adhere to the syringe nozzle.
- Gentle pressure with saline-wetted gauze can help keep the gel in place and aid clot formation.
- Additional gel can be applied if there is still evidence of haemorrhage.

Caution: Strong bleeding after haemostatic gel application may be due to misapplication. If so, follow the instructions below.

CORRECT Haemostatic Gel Application

In the event of misapplication: VETIGEL® can be removed using gauze soaked with sterile saline and/or saline lavage. Remnants of gel within the defect from the previous treatment should be fully removed, before repeating directions 3, 4 & 5 with a new haemostatic gel syringe.

64. Once haemostasis is achieved, spread the gel onto the tissue, gently pushing down and away from the centre of the gel and wound site taking care not to dislodge the gel from the wound site.

STORAGE AND HANDLING

- VETIGEL® should be stored dry and maintained between 2–25°C upon receipt.
- VETIGEL® is pre-packaged sterile and is intended for single use only.
- It is recommended that VETIGEL® be used as soon as the package is opened and unused contents discarded.

MANUFACTURER

VETIGEL® is manufactured and packaged by Cresilon, Inc., a New York City-based medical device company:

Cresilon, Inc.
86 34th Street
Suite D603/D604
Brooklyn, NY 11232
Main: +1 800-VETIGEL (+1 800-838-4435)
Email: customerservice@cresilon.com
www.cresilon.com

VETIGEL® and Cresilon® are registered trademarks of Cresilon, Inc.

DEFINITIONS OF SYMBOLS

Symbol	Definition
	Use-by date (YYYY-MM-DD)
	Consult Instructions For Use
	Caution (consult Instructions For Use)
	Manufacturer
	Date of manufacture
LOT	Batch code
STERILE	Sterile
	Do not re-use
	Do not re-sterilise
	Temperature limit
R_x	US Federal Law restricts this device to sale by or on the order of a veterinarian.
	Do not use if package is damaged

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VETIGEL® tillverkas och förpackas av det New York-baserade medicintekniska företaget Cresilon, Inc.:
VETIGEL® och Cresilon® är registrerade varumärken som tillhör Cresilon, Inc.